

INI-MOK



INJECTION-EXTRACTION FITTED CARPET CLEANER

- Suitable for injection/extraction system
- For synthetic fibres and wool
- For all types of stain

Foamless detergent for cleaning carpets and rugs specially formulated for use with an injection/extraction machine.

HOW TO USE

Prepare the solution and spray it on the surface to be cleaned. Let it act and rinse with the injection-extraction machine.



DILUTION

5% (500 ml in 10 L of water).

TECHNICAL DATA

ASPECT: Opalescent liquid

COLOR: Colorless/yellow

PERFUME: Fresh

COV: 0,30%

PH	Active substance [%]	Density (g/ml)	Pressure (bar)	Viscosity (cP)	Total dry residue [%]
9,0 ± 0,5	13 ± 1	1.025			

INI-MOK

⚠ WARNINGS

Causes serious eye irritation. Harmful to aquatic life with long lasting effects. Wash hands thoroughly after handling. Wear eye protection. If eye irritation persists: Get medical advice/attention. Contains ZINC AMMONIA CARBONATE COMPLEX, METHYLCHLOROISOTHIAZOLINONE, METHYLISOTHIAZOLINONE. May produce an allergic reaction.



WARNING

📝 NOTES

On jute backed rug or carpet check that the backing is firmly attached. Check for colour fastness in a inconspicuous area. Keep away from freezing. Before use, please always refer to the Technical Data Sheet and the Safety Data Sheet of the product. Use the PPE specified in the Safety Data Sheet. Sutter Professional cannot be considered liable for any damage due to improper product use.

Only for professional use. Safety data sheet available on request.

📦 PACKAGING

Code	Size	Quantity	Pallet
5692	Kg 5	4	9x4=36